


Using a Defibrillator

An automatic external defibrillator (AED) will measure a patient's heart rhythm to determine whether it has a 'shockable rhythm'. If the device thinks it can make a difference, it will advise a shock, shock the patient, then give CPR instructions. This will repeat in a cycle roughly every 2 minutes. If the heart does not have a shockable rhythm, it will give audible CPR instructions only.



If a shock is advised, the device may ask you to push a button, then stand clear of the patient. Follow the CPR instructions given by the device, until it warns you that it is ready to shock again.

Continue the cycle until help arrives.

ProTrainings Europe Limited

650 The Crescent, Colchester Business Park, Colchester, Essex CO4 9YQ Telephone: 01206 805359 • Email: info@protrainings.eu **www.protrainings.eu**

